Petrifilm™

6490/6491

-
-  (GB) **Staph Express Count System**
 -  (FR) **Système de numération Staph Express**
 -  (DE) **Staph Express Zählsystem**
 -  (IT) **Sistema per il conteggio Staph Express**
 -  (ES) **Sistema de recuento Staph Express**
 -  (NL) **Staph Express Count System**
 -  (SE) **Staph Express beräkningssystem**
 -  (FI) **Staph Express Count System -pikakasvatusalusta**
 -  (PT) **Sistema de Contagem Expressa de Staphylococcus**
 -  (GR) **Σύστημα Ταχείας Καταμέτρησης Staph**
 -  (JP) **黄色ブドウ球菌エクスプレス測定システム**

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AOAC Official Method 2003.07

3M Petrifilm Staph Express Plate Method for the Enumeration of *Staphylococcus aureus* in Selected Types of Processed and Prepared Foods

AOAC Official Method 2003.08

3M Petrifilm Staph Express Plate Method for the Enumeration of *Staphylococcus aureus* in Selected Dairy Foods

AOAC Official Method 2003.11

3M Petrifilm Staph Express Plate Method for the Enumeration of *Staphylococcus aureus* in Selected Meat, Seafood, and Poultry

For specific details, refer to the specific AOAC® Official MethodSM.

Nordic System for validation of alternative microbiological methods

Nordval Validation (Ref. No.: 2003-20-5408-00019) for all foods

Refer to NordVal validation for Petrifilm Staph Express Plate method details.

AFNOR validated method in compliance with ISO 16140 in comparison to ISO 6888-1 (3M – 01/9 – 04/03A) and ISO 6888-2 (3M – 01/9 – 04/03B)

Scope of the validations: All foods.

Use the following details when implementing the above Instructions for Use:

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1. FDA. 1998. Bacteriological Analytical Manual, 8th ed., (Revision A), Appendix 3.64
2. International Standards Organization, ISO 6887-1:1999, Microbiology – General guidance for the preparation of dilutions for microbiological examination
3. International Standards Organization, ISO 8261:2001, Milk and Milk Products – Preparation of test samples and dilutions for microbiological examination

AOAC is a registered trademark of AOAC INTERNATIONAL

Official Methods is a service mark of AOAC INTERNATIONAL

AFNOR is the trademark of Association française de normalization

NordVal is established as a Nordic System for validation of alternative microbiological methods (test kits).

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





9 - 15 Chilvers Road
Thornleigh, NSW 2120
Australia
1300 363 878

**3M Microbiology**

St. Paul, MN 55144-1000
Made in U.S.A. of globally sourced materials.

3M and Petrifilm are trademarks of 3M Company
© 3M 2009. All rights reserved. 34-8704-4956-7

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3M Deutschland GmbH